

Visible Implant Elastomer Manual Injection Kits

The Visible Implant Elastomer (VIE) Manual Injection Kits include all materials and supplies needed to conduct tagging programs ranging from a few individuals to thousands. For projects regularly involving many thousands of tags we recommend you consider the Air Driven Elastomer Injection System. The number of tags that can be created with each kit depends upon the size of the tags and the efficiency of usage, but as a general guide we suggest that an efficient operation might generate 200-300 tags per ml. When the supplies in your kit have been exhausted, Northwest Marine Technology offers Refill Kits also based on units of VIE.

The **60 ml Kit** contains 60 ml of VIE (customer selects up to 10 colors), two Manual Elastomer Injectors, one VI Light, one VIE Color Standard, a field carrying case, and instructions.




The **24 ml Kit** contains 24 ml of VIE (customer selects up to 8 colors), one Manual Elastomer Injector, one VI Light, one VIE Color Standard, a field carrying case, and instructions.

The **6 ml Kit** contains 6 ml of VIE (customer selects up to two colors), one Manual Elastomer Injector, one VI Light, one VIE Color Standard, a field carrying case, and instructions.



The **Trial Pack** is ideal for small projects or for pilot studies evaluating the tags for larger projects. The kit contains 1 ml of elastomer (customer selects the color), mixing and injection supplies, and instructions.



VIE Documentation

Instruction sheets, manuals and application notes are all available at www.nmt.us.

Product Sheets for VIE System, VIE Manual Injection Kits and Air Driven Elastomer Injection System.

Instruction Sheets for Manual Injection Kits

Instruction Manual for Air-Driven Elastomer Injection System

Application Notes

- Visible Implant Elastomer Material Safety Data Sheet
- Tagging small fish with Visible Implant Fluorescent Elastomer
- Instructions for tagging shrimp with the Visible Implant Elastomer (VIE) System
- Visible Implant Elastomer Project Manual

Contacting NMT

Sales and Technical Support

Northwest Marine Technology
P.O. Box 427
Shaw Island, WA 98286, USA
Telephone: (360) 468-3375
FAX: (360) 468-3844
E-mail: office@nmt.us
techsupport@nmt.us

AutoFish System

Mr. David Knutzen
955 Malin Lane SW Suite B
Tumwater, WA 98501, USA
Telephone: (360) 596-9400
FAX: (360) 596-9405
E-mail: dave.knutzen@nmt.us

Biological Support

955 Malin Lane SW Suite B
Tumwater, WA 98501, USA
Telephone: (360) 596-9400
FAX: (360) 596-9405
E-mail: biology@nmt.us

Asian Sales and Services

Mr. Yong Huang
Enfotran, Inc.
1247 106th Place NE
Bellevue, WA 98004, USA
Telephone: (425) 456-0101
FAX: (425) 456-0303
E-mail: yong.huang@enfo.us

www.nmt.us